

Apollo Silicone Cleaners | Disinfectants

Accel Five TB
Accel Prevention Concentrate
Asepticare TB-II
Birex SE
Caviwipes
Citrace Germicide
Clorox Healthcare Bleach Germicidal Spray
Clorox Disinfecting Wipes
Clorox Healthcare Bleach Germicidal Wipes
Clorox Healthcare Hydrogen Peroxide Wipes
Clorox Hydrogen Peroxide Disinfectant
Fantastik
Formula 409
Green Works All Purpose cleaner
Kics
Lysol Spray
Lysol Wipes
Oasis 146
Optim 33TB
Oxivir TB
Oxivir TB
Oxycide
PDI Sani Cloth HB
PDI Super Sani Cloth AF3
PDI Super Sani Cloth w/Bleach
Perisept
Rescue Wipes
Sani Cloth Plus
Steris Coverage Spray
Virex II 256
Virox 5 RTU
Wex-Cide 128

After cleaning, rinse area with clean water then dry with a clean, lint-free cloth. Failure to remove chemical residue will cause discoloration and damage to the material surface.